

Human GNAI1 protein (Recombinant) (His-Tag) (STJP017350)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	42.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2770
Gene Symbol	GNAI1
UniProt ID	GNAI1_HUMAN
Immunogen Region	1-354aa
Immunogen Sequence	MGCTLSA EDKAAVERSK MIDRNLRDGG EKAAREVKLL LLGAGESGKS TIVKQMKIIH EAGYSEEECK QYKAVVYSNT IQSIIAIIIRA MGRLEKIDFGD SARADDARQL FVLGAEEG FMTAELAGVI KRLWKDSGVQ ACFNRSREYQ LNDAAAYYLN DLDRIAQPNY IPTQQDVLRT RVKTTGIVET HFTFKDLHFK MFDVGG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel showing recombinant human GNAI1 protein 3µg under reducing conditions with Coomassie blue stain