

Human Cyclophilin-like 3/PPIL3 protein (Recombinant) (His-Tag) (STJP017356)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Enzyme Activity
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 0.1M NaCl, 1mM DTT
Molecular Weight	20.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	53938
Gene Symbol	PPIL3
UniProt ID	PPIL3_HUMAN
Immunogen Region	1-161aa
Immunogen Sequence	MSVTLHTDVG DIKIEVFCER TPKTCENFLA LCASNYNGC IFHRNIKGFM VQTGDPTGTG RGGNSIWGKK FEDEYSEYLG HNVRGVVSM ANGPNTNGSQ FFITYGKQPH LDMKYTVFGK VIDGLETLDE LEKLPVNEKT YRPLNDVHIK DITIHPFA Q

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Cyclophilin-like 3 PPIL3 protein 3µg reducing conditions Coomassie stain