

Human C8G protein (Recombinant) (His-Tag) (STJP017360)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 2mM DTT, 0.15M NaCl.
Molecular Weight	22.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	733
Gene Symbol	C8G
UniProt ID	CO8G_HUMAN
Immunogen Region	21-202aa
Immunogen Sequence	QKPQRPRRP ASPISTIQPK ANFDAQQFAG TWLLVAVGSA CRFLQEQGHR AEATTLHVAP QGTAMAVSTF RKLDGICWQV RQLYGDTGVL GRFLQARGA RGAVHVVVAE TDYQSFVLY LERAGQLSVK LYARSLPVSD SVLSGFQQRV QEHLTEDQI FYFPKYGFCE AADQFHLVDE VRR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human C8G protein 3µg under reducing conditions Coomassie blue stain