

Human ARF4L/ARL4D protein (Recombinant) (His-Tag) (STJP017377)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 2mM DTT, 20% Glycerol, 200mM NaCl
Molecular Weight	24.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	379
Gene Symbol	ARL4D
UniProt ID	ARL4D_HUMAN
Immunogen Region	1-201aa
Immunogen Sequence	MGNHLTEMAP TASSFLPHFQ ALHVVVIGLD SAGKTSLLYR LKFKEFVQSV PTKGFNTEKI RVPLGGSRGI TFQVWDVGGQ EKLRLPWSY TRRTDGLVFV VDAAEAERLE EAKVELHRIS RASDNQGVVPV LVLANKQDQP GALSAAEVEK RLAVRELAAL TLTHVQGCSA VDGLGLQQGL ERLYEMILKR KKAARGGKKR R

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human ARF4L ARL4D protein 3 µg Coomassie blue stain