

Human ODC1 protein (Recombinant) (His-Tag) (STJP017378)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES

Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% Glycerol, 0.1M NaCl
Molecular Weight	53.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4953
Gene Symbol	ODC1
UniProt ID	DCOR_HUMAN
Immunogen Region	1-461aa
Immunogen Sequence	MNNFGNE EFDCHFLDEG FTAKDILDQK INEVSSDDK DAFYVADLGD ILKKHLRWLK ALPRVTPFYA VKCNDKAIV KTLAATGTGF DCASKTEIQL VQSLGVPPER IYANPCKQV SQIKYAANNG VQMMTFDSEV ELMKVARAHP KAKLVLRAT DDSKAVCRLS VKFGATLRTS RLLLERAKEL NIDVVGVSFH VGSGCT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human ODC1 protein under reducing conditions Coomassie blue stain