

Human p27/Kip1 protein (Recombinant) (His-Tag) (STJP017388)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol
Molecular Weight	24.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1027
Gene Symbol	CDKN1B
UniProt ID	CDN1B_HUMAN
Immunogen Region	1-198aa
Immunogen Sequence	MSNVRVSNGS PSLERMDARQ AEHPKPSACR NLFGPVDHEE LTRDLEKHCR DMEEASQRKW NFDQNHKPL EGKYEWQVEE KGSLPEFYR PPRPPKGACK VPAQESQDVS GSRPAAPLIG APANSEDTLH VDPKTDPSDS QTGLAEQCAG IRKRPATDD S TQNKRANRT EENVSDGSPN AGSVEQTPKK PGLRRRQT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human p27 Kip1 protein under reducing conditions stained with Coomassie blue