

Human DNAJB8 protein (Recombinant) (His-Tag) (STJP017389)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 40% Glycerol, 0.15M NaCl, 1mM DTT
Molecular Weight	27.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	165721
Gene Symbol	DNAJB8
UniProt ID	DNJB8_HUMAN
Immunogen Region	1-232aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MANYEVLGV QASASPEDIK KAYRKLALRW HPDKNPDNKE EAEKKFKLVSEAYEVLSDSK KRSLYDRAGC DSWRAGGGAS TPYHSPFDTG YTFRNPEDIF REFFGGGLDPF SFEFWDSPFN SDRGGRGHGL RGAFSAGFGE FPAFMEAFSS FNMLGCSGGS HTTFSSTSFG GSSSGSSGFK SVMSSTEMIN GHKVTTKRIV ENGQERVE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human DNAJB8 protein under reducing conditions stained with Coomassie blue