

Human FKBP12.6/FKBP1B protein (Recombinant) (His-Tag) (STJP017390)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Enzyme Activity
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 1mM DTT
Molecular Weight	14.1kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2281
Gene Symbol	FKBP1B
UniProt ID	FKB1B_HUMAN
Immunogen Region	1-108aa
Immunogen Sequence	MGVEIETI SPGDGRTPFK KGQTCVVHYT GMLQNGKKFD SSRDRNKPFK FRIGKQEVK GFEEGAAQMS LGQRAKLTCT PDVAYGATGH PGVIPP NATL IFDVLLNLE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human FKBP12.6 FKBP1B protein under reducing conditions Coomassie blue stain