

Human eIF-3 alpha/EIF3J protein (Recombinant) (His-Tag) (STJP017396)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 2mM DTT, 10% Glycerol, 200mM NaCl
Molecular Weight	24 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8669
Gene Symbol	EIF3J
UniProt ID	EIF3J_HUMAN
Immunogen Region	70-258aa
Immunogen Sequence	KISEKHKIA EKIKEKERQQ KKRQEEIKR LEEPEEPKVL TPEEQADKL RLKQLQEEESD LELAKETFGV NNAVYGDAM NPSSRDDFTE FGKLLKDKIT QYEKSLYYAS FLEVLRDVC ISLEIDDLKK ITNSLTVLCS EKQKQEKQSK AKKKKKGVP GGGLKATMKD DLADYGGYDG GYVQDYEDFM

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human eIF-3 alpha EIF3J protein under reducing conditions with Coomassie blue staining