

Human Adenylosuccinate lyase/ADSL protein (Recombinant) (His-Tag) (STJP017399)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 40% Glycerol, 0.1M NaCl
Molecular Weight	59 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	158
Gene Symbol	ADSL
UniProt ID	PUR8_HUMAN
Immunogen Region	1-484aa
Immunogen Sequence	MAAG GDHGSPDSYR SPLASRYASP EMCVFVSDRY KFRTWRQLWL WLAEAEQTLG LPITDEQIQE MKSNLENIDF KMAAEKRL RHDVMAHVHT FGHCCPKAAG IHLGATSCY VGDNTDLIL RNALDLLPK LARVISRLAD FAKERASLPT LGFTHFQPAQ LTTVGKRCCL WIQDLCMDLQ NLKRVK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human ADSL protein 3µg reducing conditions Coomassie blue stain