

Human GADD45GIP1 protein (Recombinant) (His-Tag) (STJP017404)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 40% Glycerol, 2mM DTT, 0.2M NaCl.
Molecular Weight	22.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	90480
Gene Symbol	GADD45GIP1
UniProt ID	G45IP_HUMAN
Immunogen Region	48-222aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MPRWQLGPRY AAKQFARYGA ASGVVPGSLW PSPEQLRELE AEEREWYPSL ATMQESLRVK QLAAEQKRRE REQHIAECMA KMPQMIWNWQ QQQRENWEKA QADKERRARL QAFAQELLGY QVDPRSARFQ ELLQDLEKKE RKRLKEEKQK RKKEARAAAL AAFAQDPAA SGAPSS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human GADD45GIP1 protein under reducing conditions stained with Coomassie blue