

Human GMDS protein (Recombinant) (His-Tag) (STJP017406)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 30% Glycerol, 0.1M NaCl, 1mM DTT, 0.1mM PMSF
Molecular Weight	44.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2762
Gene Symbol	GMDS
UniProt ID	GMDS_HUMAN
Immunogen Region	1-372aa
Immunogen Sequence	MAHAPARCPs ARGSGDGEMG KPRNVALITG ITGQDGsYLA EFLEKGYEV HGIVRRSSsF NTGRIEHLyK NPQAHIEGNM KLHYGDLTDS TCLVKIINEV KPTEIYNLGA QSHVKISFDL AEYTADVDGV GTLRLLDAVK TCGLINSVKF YQASTSELYG KVQEIPQKET TPFYPRSPYG AAKLYAYWIV VNFREAYNLF AVNGIL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human GMDS protein stained with Coomassie blue