

Human I-TAC/CXCL11 protein (Recombinant) (His-Tag) (STJP017420)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 10mM Sodium Citrate buffer (pH 3.5) containing 2mM DTT, 20% Glycerol
Molecular Weight	10.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6373
Gene Symbol	CXCL11
UniProt ID	CXL11_HUMAN
Immunogen Region	22-94aa
Immunogen Sequence	FPMFKRGRC LCIGPGVKAV KVADIEKASI MYPSNNCDKI EVIITLKENK GQRCLNPKSK QARLIKKVE RKNF

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human I-TAC CXCL11 protein under reducing conditions with Coomassie blue stain