

Human p38 beta/MAPK11 protein (Recombinant) (His-Tag) (STJP017423)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 2mM DTT, 20% Glycerol, 100mM NaCl
Molecular Weight	43.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5600
Gene Symbol	MAPK11
UniProt ID	MK11_HUMAN
Immunogen Region	1-364aa
Immunogen Sequence	MSGPRAG FYRQELNKTV WEVPQRLQGL RPYGSGAYGS VCSAYDARLR QKVAVKKLSR PFQSLIHARR TYRELRLKKH LKHENVIGLL DVFTPATSIIE DFSEVYLVT LMGADLNNIV KCQALSDEHV QFLVYQLLRG LKYIHSAGII HRDLKPSNVA VNEDCELRIL DFGLARQADE EMTGYVATRW YRAPEIMLNW MHYNQT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human p38 beta MAPK11 protein under reducing conditions stained with Coomassie blue