

Human PSMB4 protein (Recombinant) (His-Tag) (STJP017430)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 30% Glycerol, 0.1M NaCl
Molecular Weight	26.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5692
Gene Symbol	PSMB4
UniProt ID	PSB4_HUMAN
Immunogen Region	46-264aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MTQNPMTGT SVLGVKFEGG VVIAADMLGS YGSLARFRNI SRIMRVNNST MLGASGDYAD FQYLKQVLGQ MVIDEELLGD GHSYSPRAIH SWLTRAMYSR RSKMNPWNT MVIIGGYADGE SFLGYVDMLG VAYEAPSLAT GYGAYLAQPL LREVLEKQPV LSQTEARDLV ERCMRVLYYR DARSYNRFQI ATVTEKGVEI EGPLSTET

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human PSMB4 protein under reducing conditions with Coomassie blue stain