

Human Importin beta/KPNB1 protein (Recombinant) (His-Tag) (STJP017440)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 30% Glycerol, 0.1M NaCl
Molecular Weight	99.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3837
Gene Symbol	KPNB1
UniProt ID	IMB1_HUMAN
Immunogen Region	1-876aa
Immunogen Sequence	MELITIL EKTVSPDRLE LEAAQKFLER AAVENLPTFL VELSRVLNAP GNSQVARVAA GLQIKNSLTS KDPDIKAQYQ QRWLAI DANA RREVKNYVLQ TLGTETYRPS SASQCVAGIA CAEIPVNQWP ELIPQLVANV TNPNSTEHMK ESTLEAIGYI CQDIDPEQLQ DKSNEILTAI IQGMRKEEPS NNVKLAATNA LLNSLE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Importin beta KPNB1 protein under reducing conditions