

Human PSMA8 protein (Recombinant) (His-Tag) (STJP017444)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% Glycerol, 0.1M NaCl
Molecular Weight	30 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	143471
Gene Symbol	PSMA8
UniProt ID	PSMA8_HUMAN
Immunogen Region	1-250aa
Immunogen Sequence	MGSSHHHHH SSGLVPRGSH MASRYDRAIT VFSPDGHLEFQ VEYAEAVKK GSTAVGIRGT NIVVLGVEKK SVAKLQDERT VRKICALDDH VCMAFAGLTA DARVVINRAR VECQSHKLTVE EDPVTVEYIT RFIATLKQKY TQSNRRPFG ISALIVGFDD DGISRLYQTD PSGTYHAWKA NAIGRSKTV REFLEKNYTE DAIASDSEAI KLAIKALLEV VQSGGKNI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel showing recombinant human PSMA8 protein 3µg under reducing conditions with Coomassie blue staining