

Human Calcium Activated Nucleotidase 1/CANT1 protein (Recombinant) (His-Tag) (STJP017461)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% Glycerol, 50mM NaCl
Molecular Weight	40.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	124583
Gene Symbol	CANT1
UniProt ID	CANT1_HUMAN
Immunogen Region	63-401aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSHMRPAPG RPPTHNAHNW RLQAPANWY NDTYPLSPPQ RTPAGIRYRI AVIADLTES RAQEENTWFS YLKKGYLTL SSGDKVAV EW DKDHGVLESH LAEKGRGMEL SDLIVFNGKL YSVDDRTGVV YQIEGSKAVP WVILSDGDGT VEKGFKAEWL AVKDERLYVG GLGKEWTTTT GDVVNENPEW VKVVG YKGSV DHENWVSN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human CANT1 protein 3µg under reducing conditions Coomassie stain