

Human LCMT1 protein (Recombinant) (His-Tag) (STJP017469)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 1mM DTT, 0.15M NaCl.
Molecular Weight	40.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51451
Gene Symbol	LCMT1
UniProt ID	LCMT1_HUMAN
Immunogen Region	1-334aa
Immunogen Sequence	MATRQR ESSITSCCST SSCDADDEGV RGTCEDASLC KRFVSIQYW HDPYIQHFVR LSKERKAPEI NRGYFARVHG VSQLIKAFRLR KTECHCQIVN LGAGMDTTFW RLKDEDLLPS KYFEVDFPMI VTRKLHSIKC KPPLSSPILE LHSEDTLQMD GHILDSKRYA VIGADLRDLS ELEEKLKKN MNTQLPTLLI AECVLV

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant human LCMT1 protein SDS-PAGE gel 3µg reducing conditions Coomassie blue stain