

Human RTN4IP1 protein (Recombinant) (His-Tag) (STJP017483)

GENERAL INFORMATION	
Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli
PRODUCT PROPERTIES	
Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 30% Glycerol, 0.1M NaCl
Molecular Weight	41.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	84816
Gene Symbol	RTN4IP1
UniProt ID	RT4I1_HUMAN
Immunogen Region	41-396aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSTVMPAWV IDKYGKNEVL RFTQNMMMPI IHYPNEVIVK VHAASVNPID VNMRSYGAT ALNMKRDPLH VKIKGEFPL TLGRDVSGVV MECGLDVKYF KPGDEVWAAV PPWKQGTLSF FVVVSGNEVS HKPKSLTHTQ AASLPYVALT AWSAINKVG G LNDKNCTGKR VLILGASGGV GTFAIQVMKA WDAHVTAVCS QDASELVR
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RTN4IP1 protein under reducing conditions with Coomassie blue stain