

Human YEATS4 protein (Recombinant) (His-Tag) (STJP017488)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 50% Glycerol, 2mM DTT, 0.2M NaCl, 0.1mM PMSF, 50mM Imidazole
Molecular Weight	28.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8089
Gene Symbol	YEATS4
UniProt ID	YEATS4_HUMAN
Immunogen Region	1-227aa
Immunogen Sequence	MFKRMAE FGPDSGGRVK GVTIVKPIVY GNVARYFGKK REEDGHTHQW TVYVKPYRNE DMSAYVKKIQ FKLHESYGNP LRVVTKPPYE ITETGWGEFE IIKIFFIDP NERPVTLYHL LKLFQSDTNA MLGKKTIVSE FYDEMIFQDP TAMMQQLTT SRQLTLGAYK HETEF AELEV KTRKLEAAK KKTSFEIAEL KERLKA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human YEATS4 protein stained with Coomassie blue