

Human RACK1 protein (Recombinant) (His-Tag) (STJP017498)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 50% Glycerol, 0.2M NaCl, 5mM DTT
Molecular Weight	37.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10399
Gene Symbol	RACK1
UniProt ID	RACK1_HUMAN
Immunogen Region	1-317aa
Immunogen Sequence	MTEQMTLRGT LKGHNGWVTQ IATTPQFPDM ILSASRDKTI IMWKLTRDET NYGIPQRALR GHSHFVSDVV ISSDGQFALS GSWDGTLLRLW DLTTGTTTTRR FVGHTKDVLS VAFSSDNRQI VSGSRDKTIK LWNTLGVCCKY TVQDESHSEW VSCVRFSPNS SNPIIVSCGW DKLVKVNLA NCKLKTNHIG HTGYLNTVTV SPDGSL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RACK1 protein 3 µg under reducing conditions with Coomassie stain