

Human CYB5R3 protein (Recombinant) (His-Tag) (STJP017503)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 1mM DTT, 0.1M NaCl.
Molecular Weight	34 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1727
Gene Symbol	CYB5R3
UniProt ID	NB5R3_HUMAN
Immunogen Region	27-301aa
Immunogen Sequence	FQRST PAITLESPDI KYPLRLIDRE IISHDTRFRF FALPSPQHIL GLPVGQHIYL SARIDGNLVV RPYTPISDD DKGFVDLVK VYFKDTHPKF PAGGKMSQYL ESMQIGDTIE FRGFSGLLVY QGKGF AIRP DKKSNPIRT VKSVGMIAGG TGITPMLQVI RAIMKDPDDH TVCHLLFANQ TEKDILLRPE LEELRN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human CYB5R3 protein 3µg under reducing conditions stained with Coomassie blue