

Human FGF-10 protein (Recombinant) (His-Tag) (STJP017510)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 50% Glycerol, 0.2M NaCl, 2mM DTT, 2mM EDTA
Molecular Weight	22.0 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2255
Gene Symbol	FGF10
UniProt ID	FGF10_HUMAN
Immunogen Region	38-208aa
Immunogen Sequence	QALGQ DMVSPEATNS SSSSFSSPSS AGRHVRSYNH LQGDVRWRKL FSFTKYFLKI EKNGKVS GTK KENCPSILE ITSVEIGVVA VKAINSYYL AMNKKGKLYG SKEFNNDCKL KERIEENGYN TYASFNWQHN GRQMYVALNG KGAPRRGQKT RRKNTSAHFL PMVVHS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human FGF-10 protein under reducing conditions Coomassie blue stain