

Human PAPSS1 protein (Recombinant) (His-Tag) (STJP017513)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 0.1M NaCl
Molecular Weight	70.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9061
Gene Symbol	PAPSS1
UniProt ID	PAPS1_HUMAN
Immunogen Region	24-624aa
Immunogen Sequence	RATNV TYQAHVSRN KRGQVVGTRG GFRGCTVWLT GLSGAGKTTV SMALEEYLVC HGIPCYTLDG DNIRQGLNKN LGFSPEDREE NVRRIAEVAK LFADAGLVC I TSFISPYTQD RNNARQIHEG ASLPFFFEV DAPLHVCEQR DVKGLYKKAR AGEIKGFTGI DSEYEKPEAP ELVLKTDSCD VNDCVQQVVE LLQERD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human PAPSS1 protein under reducing conditions with Coomassie blue stain