

Human RPLP2 protein (Recombinant) (His-Tag) (STJP017528)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% Glycerol, 0.1M NaCl
Molecular Weight	14.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6181
Gene Symbol	RPLP2
UniProt ID	RLA2_HUMAN
Immunogen Region	1-115aa
Immunogen Sequence	MRYVAS YLLAALGGNS SPSAKDIKKI LDSVGIEADD DRLNKVISEL NGKNIEDVIA QGIGKLASVP AGGAVAVSAA PGSAAPAAGS APAAAEKKD EKKEESEESD DDMGFGFLFD

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human RPLP2 protein 3µg under reducing conditions Coomassie blue stain