

Human GMPR protein (Recombinant) (His-Tag) (STJP017533)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 40% Glycerol, 0.15M NaCl, 1mM DTT
Molecular Weight	39.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2766
Gene Symbol	GMPR
UniProt ID	GMPR1_HUMAN
Immunogen Region	1-345aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MPRIDADLKL DFKDVLLRPK RSSLKSRAEV DLERTFTFRN SKQTYSGIPI IVANMDTVGT FEMAAMVSQH SMFTAIHKHY SLDDWKL FAT NHPECLQNVV VSSGSGQNDL EKMTSILEAV PQVKFICLDV ANGYSEHFVE FVKLVRAKFP EHTIMAGNVV TGEMVEELIL SGADIHKVGV GPGSVCTTRT KTGVGYPQLS AVIECADS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human GMPR protein under reducing conditions stained with Coomassie blue