

Human HNRNPA1 protein (Recombinant) (His-Tag) (STJP017534)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 40% Glycerol, 0.15M NaCl, 1mM DTT
Molecular Weight	36.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3178
Gene Symbol	HNRNPA1
UniProt ID	ROA1_HUMAN
Immunogen Region	1-320aa
Immunogen Sequence	MSKSESP KEPEQLRKLF IGGLSFETTD ESLRSHFEQW GTLTDCVVMR DPNTKRSRGF GFVITYATVEE VDAAMNARPH KVDGRVVEPK RAVSREDSQR PGAHLTVKKI FVGGIKEDTE EHHLRDYFEQ YGKIEVIEIM TDRGSGKKRG FAFVTFDDHD SVDKIVIQKY HTVNGHNCEV RKALSKQEMA SASSSQRGRS GSGNFG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human HNRNPA1 protein 3µg under reducing conditions Coomassie stain