

Human IκB-beta protein (Recombinant) (His-Tag) (STJP017551)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 1mM DTT
Molecular Weight	40.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4793
Gene Symbol	NFKB1B
UniProt ID	IKBB_HUMAN
Immunogen Region	1-356aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSHMAGVAC LGKAADAEW CDSGLGSLGP DAAAPGGPGL GAELGPGLSW APLVFGYVTE DGD TALHLAV IHQHEPFLDF LLGFSAGTEY MDLQNDLGQT ALHLAAILGE TSTVEKLYAA GAGLCVAERR GHTALHLACR VGAHACARAL LQPRPRRPRE APDTYLAQGP DRTPDNTNHTP VALYPDSLE KEEEESEEDW KLQLEAEN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human IκB-beta protein under reducing conditions