

Human RPL8 protein (Recombinant) (His-Tag) (STJP017565)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 2mM DTT, 50% Glycerol, 300mM NaCl, 2mM EDTA
Molecular Weight	30.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6132
Gene Symbol	RPL8
UniProt ID	RL8_HUMAN
Immunogen Region	1-257aa
Immunogen Sequence	MGRVIRGQRK GAGSVFRAHV KHRKGAARLR AVDFEAERHGY IKGIVKDIH DPGRGAPLAK VVFRDPYRFK KRTELFIAAE GIHTGQFVYC GKKAQLNIGN VLPVGTMP EG TIVCCLEEKP GDRGKLARAS GNYATVISHN PETKKTRVKL PSGSKKVISS ANRAVVG VVA GGGRIDKPIL KAGRAYHKYK AKRNCWPRVR GVAMNP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel showing recombinant human RPL8 protein 3µg under reducing conditions
Coomassie blue stain