

Human VEGF 165 protein (Recombinant) (No-Tag) (STJP017572)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH8.5) containing 2mM EDTA, 0.1mM PMSF, 1mM DTT, 10% Glycerol
Molecular Weight	19.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7422
Gene Symbol	VEGFA
UniProt ID	VEGFA_HUMAN
Immunogen Region	207-371aa
Immunogen Sequence	MAPMAEGGGQ NHHEVVKFMD VYQRSYCHPI ETLVDIFQEY PDEIEYIFKP SCVPLMRCGG CCNDEGLECV PTEESNITMQ IMRIKPHQGQ HIGEMSFLQH NKCECRPKKD RARQENPCGP CSERRKHLFV QDPQTKCSC KNTDSRCKAR QLELNERTCR CDKPRR

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human VEGF 165 protein under reducing conditions with Coomassie blue stain