

Human Casein kinase 2 alpha 1/CSNK2A1 protein (Recombinant) (His-Tag) (STJP017578)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 0.1M NaCl, 1mM DTT, 0.1mM PMSF
Molecular Weight	47.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1457
Gene Symbol	CSNK2A1
UniProt ID	CSK21_HUMAN
Immunogen Region	1-391aa
Immunogen Sequence	MSGPV PSRARVYTDV NTHRPREYWD YESHVVEWGN QDDYQLVRKL GRGKYSEVFE AINITNNEKV VVKILKPVKK KKKIKREIKIL ENLRGGPNII TLADIVKDPV SRTPALVFEH VNNTDFKQLY QTLTDYDIRF YMYEILKALD YCHSMGIMHR DVKPHNVIMID HEHRKRLRID WGLAEFYHPG QEYNVRVASR YFKGPE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human Casein kinase 2 alpha 1 CSNK2A1 protein under reducing conditions