

Human LSM3 protein (Recombinant) (His-Tag) (STJP017582)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 2mM DTT, 20% Glycerol, 100mM NaCl
Molecular Weight	14.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27258
Gene Symbol	LSM3
UniProt ID	LSM3_HUMAN
Immunogen Region	1-102aa
Immunogen Sequence	MADDVDQ QTTNTVEEP LDLIRLSLDE RIYVKMRNDR ELRGRRLHAYD QHLNMLGDV EETVTTIEID EETYEIYKS TKRNIPMLFV RGDGVVLVAP PLRVG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human LSM3 protein stained with Coomassie blue