

Human ADH5 protein (Recombinant) (His-Tag) (STJP017594)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 20% Glycerol 0.1M NaCl, 1mM DTT
Molecular Weight	42.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	128
Gene Symbol	ADH5
UniProt ID	ADHX_HUMAN
Immunogen Region	1-374aa
Immunogen Sequence	MANEVI KCKAAVAWEA GKPLSIEEIE VAPPKAHEVR IKIITAVCH TDAYTLSGAD PEGCFPVILG HEGAGIVESV GEGVTKLKAG DTVIPLYIPQ CGECKFCLNP KTNLCQKIRV TQKGGLMPDG TSRFTCKGKT ILHYMGTSSTF SEYTVVADIS VAKIDPLAPL DKVCLLGCGI STGYGAAVNT AKLEPGSVCA VFGLGG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human ADH5 protein 3 µg under reducing conditions
Coomassie blue stain