

Human LSM2 protein (Recombinant) (His-Tag) (STJP017599)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 5mM DTT, 40% Glycerol, 200mM NaCl
Molecular Weight	13.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57819
Gene Symbol	LSM2
UniProt ID	LSM2_HUMAN
Immunogen Region	1-95aa
Immunogen Sequence	MGSSHHHHHH SSSLVPRGSH MGSFMLFYF FKSIVGKDVV VELKNDLSIC GTLHSVDQYL NIKLTDISVT DPEKYPHMLS VKNCFIRGSV VRYVQLPADE VDTQLLQDAA RKEALQQKQ

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human LSM2 protein under reducing conditions stained with Coomassie blue