

Human EYA2 protein (Recombinant) (His-Tag) (STJP017601)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 40% Glycerol, 0.15M NaCl, 1mM DTT
Molecular Weight	33.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2139
Gene Symbol	EYA2
UniProt ID	EYA2_HUMAN
Immunogen Region	244-514aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMERVFVW DLDETHIIFH SLLTGTFFASR YGKDTTTSVR IGLMMEEMIF NLADTHLFFN DLEDQDIHV DDVSSDDNGQ DLSTYNFSAD GFHSSAPAAAN LCLGSGVHGG VDWMRKLAFR YRRVKEMYNT YKNNVGLIG TPKRETWLQL RAELEALTDL WLTHSLKALN LINSRPNVCVN VLVTTTQLIP ALAKVLLYGL GSVFPIEN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human EYA2 protein under reducing conditions stained with Coomassie blue