

Human ATOH1 protein (Recombinant) (No-Tag) (STJP017604)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 0.4M urea
Molecular Weight	21.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	474
Gene Symbol	ATOH1
UniProt ID	ATOH1_HUMAN
Immunogen Region	158-354aa
Immunogen Sequence	MKQRRRLAANA RERRRMHGLN HAFDQLRNVI PSFNNDKKLS KYETLQMAQI YINALSELLQ TPGGGEQPPP PPASCKSDHH HLRTAASYEG GAGNATAAGA QQASGGGSRP TPPGSCRTRF SAPASAGGYS VQLDALHFST FEDSALTAMM AQKNLSPSLP GSILQPVQEE NSKTSPRSHR SDGEFSPHSY YSDSDEAS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human ATOH1 protein under reducing conditions stained with Coomassie blue