

Human Glutathione S-transferase omega 2/GSTO2 protein (Recombinant) (His-Tag) (STJP017629)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.5) containing 0.2M NaCl, 1mM DTT, 40% Glycerol
Molecular Weight	30.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	119391
Gene Symbol	GSTO2
UniProt ID	GSTO2_HUMAN
Immunogen Region	1-243aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMGSDATR TLGKGSQPPG PVPEGLIRIY SMRFCPYSHR TRLVLKAKDI RHEVVNINLR NKPEWYYTKH PFGHIPVLET SQCLYIESV IACEYLDDAY PGRKLFYDP YERARQKMLL ELFCVKPHLT KECLVALRCG RECTNLKAAL RQEFSNLEEI LEYQNTTFFG GTCISMIDYL LWPWFERLDV YGILDCVSHT PALRLWIS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human GSTO2 protein under reducing conditions Coomassie stain