

Human BCAT2 protein (Recombinant) (His-Tag) (STJP017636)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 10% Glycerol 0.2M NaCl, 1mM DTT
Molecular Weight	43.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	587
Gene Symbol	BCAT2
UniProt ID	BCAT2_HUMAN
Immunogen Region	28-392aa
Immunogen Sequence	ASSSF KAADLQLEMT QKPHKKPGPG EPLVFGKTFT DHMLMVEWND KGWGQPRIQP FQNLTLPAS SSLHYSLQLF EGMKAFKGD QQVRLFRPWL NMDRMLRSAM RLCLPSFDKL ELLECIRRLI EVDKDWVPDA AGTSLYVRPV LIGNEPSLGV SQPTRALLFV ILCPVGAYFP GGSVTPVSLI ADPAFIRAWV GGVGNY

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human BCAT2 protein stained with Coomassie blue