

Human GTF2F2 protein (Recombinant) (His-Tag) (STJP017637)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol 0.2M NaCl, 1mM DTT
Molecular Weight	30.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2963
Gene Symbol	GTF2F2
UniProt ID	T2FB_HUMAN
Immunogen Region	1-249aa
Immunogen Sequence	MAERGELDLT GAKQNTGVWL VKVPKYLSSQ WAKASGRGEV GKLRIAKTQG RTEVSFTLNE DLANIHDIGG KPASVSAPRE HPFVLQSVGG QTLTVFTESS SDKLSLEGIV VQRAECRPAA SENYMLKRL QIEESSKPVR LSQQLDKVVT TNYKPVANHQ YNIEYERKKK EDGKRARADK QHVLDMFLSA FEKHQYYNLK DLVDIT

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human GTF2F2 protein under reducing conditions Coomassie blue stain