

Human NEIL2 protein (Recombinant) (His-Tag) (STJP017639)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% Glycerol, 0.1M NaCl
Molecular Weight	39.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	252969
Gene Symbol	NEIL2
UniProt ID	NEIL2_HUMAN
Immunogen Region	1-332aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSHMPEGPL VRKFHHLVSP FVGQQVVKTG GSSKKLQPAS LQSLWLQDTQ VHGGKFLRF DLDEEMGPPG SSPTPEPPQK EVQKEGAADP KQVGEPGSGQK TLDGSSRS AE LVPGEDDSE YLERDAPAGD AGRWLRVSFG LFGSVVWVND F SRAKKANKRG DWRDPSPRLV LHFGGGGFLA FYNCQLSWSS SPVVTPTCDI LSEKFHRG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel showing recombinant human NEIL2 protein 3ug under reducing conditions with Coomassie blue stain