

## Human HOMER3 protein (Recombinant) (His-Tag) (STJP017659)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Applications  | SDS-PAGE |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                    |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| Concentration       | 0.25 mg/mL                                                                                                                                         |
| Formulation         | Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 30% Glycerol, 0.15M NaCl, 1mM DTT                                                               |
| Molecular Weight    | 42.5 kDa                                                                                                                                           |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                                       |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">9454</a>                                                                                                                                                                                                                  |
| Gene Symbol        | <a href="#">HOMER3</a>                                                                                                                                                                                                                |
| UniProt ID         | <a href="#">HOME3_HUMAN</a>                                                                                                                                                                                                           |
| Immunogen Region   | 1-361aa                                                                                                                                                                                                                               |
| Immunogen Sequence | MSTAR EQPIFSTRAH VFQIDPATKR NWIPAGKHAL TVSYFYDATR NVYRIISIGG AKAIINSTVT<br>PNMTFTKTSQ KFGQWADSR A NTVYGLGFAS EQHLTQFAEK FQEVKEAARL AREKSQDGG E LTSPALGLAS<br>HQVPPSP LVS ANGPGEELF RSQSADAPGP TERERLKKML SEGSGVEVQW EAEFFALQDS NNKLAG |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant human HOMER3 protein SDS-PAGE gel 3 µg Coomassie blue stain reducing conditions