

Human JAM-A/F11R protein (Recombinant) (His-Tag) (STJP017666)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 0.15M NaCl, 1mM DTT
Molecular Weight	25.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	50848
Gene Symbol	F11R
UniProt ID	JAM1_HUMAN
Immunogen Region	26-238aa
Immunogen Sequence	LGSVT VHSSEPEVRI PENNPVKLSC AYSGFSSPRV EWKFDQGD TT RLVCYNNKIT ASYEDRVTF LPTGITFKSVT REDTGTYTCM VSEEGGNSYG EVKVKLIVLV PPSKPTVNIP SSATIGNRAV LTCSEQDGSP PSEYTWFKDG IVMPNPKST RAFSNSSYVL NPPTTGELVFD PLSASDTGEY SCEARNGYGT PMTSNA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human JAM-A F11R protein under reducing conditions with Coomassie stain