

Human CLTA protein (Recombinant) (His-Tag) (STJP017670)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% Glycerol
Molecular Weight	26.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1211
Gene Symbol	CLTA
UniProt ID	CLCA_HUMAN
Immunogen Region	1-218aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSMAELDP FGAPAGAPGG PALGNGVAGA GEEDPAAAF AQQESEIAGI ENDEAFALD GGAPGPQPHG EPPGGPDAVD GVMNGEYYQE SNGPTDSYAA ISQVDRQLQE PESIRKWREE QMERLEALDA NSRKQAEWK EKAIKELEEW YARQDEQLQK TKANNRAAEE AFVNDIDESS PGTEWERVER LCDFNPKSSK QAKDVSRM

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human CLTA protein under reducing conditions Coomassie blue stain