

Human Kallikrein 3/PSA protein (Recombinant) (His-Tag) (STJP017676)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Denatured
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.4M Urea, 10% Glycerol
Molecular Weight	28.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	354
Gene Symbol	KLK3
UniProt ID	KLK3_HUMAN
Immunogen Region	25-261aa
Immunogen Sequence	IVGGW ECEKHSQPWQ VLVASRGRAV CGGVLVHPQW VLTAAHCI RN KSVILLGRHS LFHPEDTGQV FQVSHSFPH LYDMSLLKNR FLRPGDDSSH DLMMLRLSEP AELTDAVKVM DLPTQEPALG TTCYASGWGS IEPEEF LTPK KLQCVDLHVI SNDVCAQVHP QKVTKFMLCA GRWTGGKSTC SGDSGGPLVC NGVLQG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Kallikrein 3
PSA protein stained with Coomassie blue