

Human MSRA protein (Recombinant) (His-Tag) (STJP017679)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 1mM DTT, 50mM NaCl
Molecular Weight	26.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4482
Gene Symbol	MSRA
UniProt ID	MSRA_HUMAN
Immunogen Region	24-235aa
Immunogen Sequence	MGSSHHHHH SGLVPRGSH MGSHEMGNAS NIVSPQEALP GRKEQTPVAA KHHVNGNRTV EPFPEGTQMA VFGMGCFWGA ERKFWVLKGV YSTQVGFAGG YTSNPTYKEV CSEKTGHAEV VRVYQPEHM SFEELLKVFV ENHDPTQGMR QGNDHGTQYR SAIYPTSAKQ MEALSSKEN YQVLSEHGF GPITTDIREG QTFYYAEDYH QQYLSKNPNG YCGLGGTG

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human MSRA protein under reducing conditions stained with Coomassie blue