

## Human NRK1/NMRK1 protein (Recombinant) (His-Tag) (STJP017680)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Applications  | SDS-PAGE |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                    |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| Concentration       | 0.5 mg/mL                                                                                                                                          |
| Formulation         | Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol                                                                               |
| Molecular Weight    | 25.6 kDa                                                                                                                                           |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                              |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">54981</a>                                                                                                                                                                                                        |
| Gene Symbol        | <a href="#">NMRK1</a>                                                                                                                                                                                                        |
| UniProt ID         | <a href="#">NRK1_HUMAN</a>                                                                                                                                                                                                   |
| Immunogen Region   | 1-199aa                                                                                                                                                                                                                      |
| Immunogen Sequence | MKTFIIG ISGVTNSGKT TLAKNLQKHL PNCVISQDD FFKPESEIET DKNGLQYDV LEALNMEKMM<br>SAISCWMESA RHSVSTQDE SAEIIPILII EGFLLFNYKP LDTIWNRSYF LTIPYEECKR RRSTRVYQPP<br>DSPGYFDGHV WPMYLYRQE MQDITWEVVY LDGKSEEDL FLQVYEDLIQ ELAKQKCLQV TA |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human NRK1 NMRK1 protein 3ug under reducing conditions