

Human VCAM-1/CD106 protein (Recombinant) (His-Tag) (STJP017683)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	76 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7412
Gene Symbol	VCAM1
UniProt ID	VCAM1_HUMAN
Immunogen Region	25-698aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MFKIETTPES RYLAQIGDSV SLTCSTTGCE SPFFSWRTQI DSPLNGKVTN EGTTSTLTMN PVSFNGNEHSY LCTATCESRK LEKGIQVEIY SFPKDPPEIHL SGPLEAGKPI TVKCSVADVY PFDRLEIDLK KGDHLMKSQE FLEDADRKSL ETKSLEVTFT PVIEDIGKVL VCRACLHIDE MDSVPTVRQA VKELQVYISP KNTVISVN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human VCAM-1 CD106 protein under reducing conditions