

Human TBCB protein (Recombinant) (His-Tag) (STJP017691)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 20% Glycerol
Molecular Weight	29.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1155
Gene Symbol	TBCB
UniProt ID	TBCB_HUMAN
Immunogen Region	1-244aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MGSHEVTGV SAPTVTVFIS SSLNTRFRSEK RYSRSLTIAE FKCKLELLVG SPASCMEL EL YGVDDKFYSK LDQEDALLGS YPVDDGCRH VIDHSGARLG EYEDVSRVEK YTISQEAYDQ RQDTVRSFLK RSKLGRYNEE ERAQQEAEAA QRLAEEKAQA SSIPVGSRC E VRAAGQSPRR GTVMYVGLTD FKPGYWIGVR YDEPLGKN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human TBCB protein under reducing conditions with Coomassie blue stain