

Human Cyclin B2 protein (Recombinant) (His-Tag) (STJP017693)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 10% Glycerol, 0.5M NaCl, 250mM Imidazol
Molecular Weight	47.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9133
Gene Symbol	CCNB2
UniProt ID	CCNB2_HUMAN
Immunogen Region	1-398aa
Immunogen Sequence	HMALLRR PTVSSDLENI DTGVNSKVKS HVTIRRTVLE EIGNRVTTTRA AQVAKKAQNT KVPVQPTKTT NVNKQLKPPTA SVKPVQMEKL APKGPSPTPE DVSMKEENLC QAFSDALLCK IEDIDNEDWE NPQLCSDYVK DIYQYLRQLE VLQSPNPHFL DGRDINGRMR AILVDWLQVQV HSKFRLLQET LYMCVGIMDR FLQVQP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Cyclin B2 protein under reducing conditions stained with Coomassie blue